


PERIPHERAL

HeliFlex^{Ti}

1ST PERIPHERAL SELF-EXPANDING

Bio **ACTIVE** STENT




HEXACATH

HeliFlex^{Ti}

1ST PERIPHERAL SELF-EXPANDING Bio ACTIVE STENT




LOWER RESTENOSIS



Bio Active Coating (TiNO)

- Reduces platelet and fibrin deposition¹
- Minimizes corrosion and inflammation²
- Promotes re-endothelialization³
- Decreases thrombus formation⁴

SUPERIOR FRACTURE RESISTANCE



Helicoidal Design

- Virtually fracture free
- Unmatched flexibility
- Highest conformability



ENHANCED VISIBILITY



Tantalum Markers

- Optimal stent visibility
- Accurate stent positioning

EXTENSIVE RANGE



**EASY TO USE
PULL BACK SYSTEM**

PIONEER IN BIO ACTIVE COATING

HEXACATH

Heliflex^{Ti}

PERIPHERAL SELF-EXPANDING BIO ACTIVE STENT (B.A.S)

- Usable lengths:
80cm or 140cm
- Minimum introducer sheath:
5F: for Ø4 & 5mm
6F: for Ø6 to 10mm
- Stent integrated tantalum markers
at each extremity:
3 for Ø4 & 5mm
5 for Ø6 to 10mm
- Stent material composition:
Nitinol
- Coating: Titanox™
(Titanium-Nitride-Oxide)

OTW	OTW
0.035"	0.018"

Stent Length (mm)	Stent Diameter (mm)	OTW 0.018"	OTW 0.035"	Catheter Lengths	
				80cm Reference	140cm Reference
30	4	●		HFXI 4-30 Ti L80	HFXI 4-30 Ti L140
	5	●		HFXI 5-30 Ti L80	HFXI 5-30 Ti L140
	6		●	HFXI 6-30 Ti L80	HFXI 6-30 Ti L140
	7	●		HFXI 7-30 Ti L80	HFXI 7-30 Ti L140
	8	●		HFXI 8-30 Ti L80	HFXI 8-30 Ti L140
	9	●		HFXI 9-30 Ti L80	HFXI 9-30 Ti L140
40	4	●		HFXI 4-40 Ti L80	HFXI 4-40 Ti L140
	5	●		HFXI 5-40 Ti L80	HFXI 5-40 Ti L140
	6		●	HFXI 6-40 Ti L80	HFXI 6-40 Ti L140
	7	●		HFXI 7-40 Ti L80	HFXI 7-40 Ti L140
	8	●		HFXI 8-40 Ti L80	HFXI 8-40 Ti L140
	9	●		HFXI 9-40 Ti L80	HFXI 9-40 Ti L140
60	4	●		HFXI 4-60 Ti L80	HFXI 4-60 Ti L140
	5	●		HFXI 5-60 Ti L80	HFXI 5-60 Ti L140
	6		●	HFXI 6-60 Ti L80	HFXI 6-60 Ti L140
	7	●		HFXI 7-60 Ti L80	HFXI 7-60 Ti L140
	8	●		HFXI 8-60 Ti L80	HFXI 8-60 Ti L140
	9	●		HFXI 9-60 Ti L80	HFXI 9-60 Ti L140
80	6		●	HFXI 6-80 Ti L80	HFXI 6-80 Ti L140
	7	●		HFXI 7-80 Ti L80	HFXI 7-80 Ti L140
	8	●		HFXI 8-80 Ti L80	HFXI 8-80 Ti L140
	9	●		HFXI 9-80 Ti L80	HFXI 9-80 Ti L140
	10	●		HFXI 10-80 Ti L80	HFXI 10-80 Ti L140
	100	6		●	HFXI 6-100 Ti L80
7		●		HFXI 7-100 Ti L80	HFXI 7-100 Ti L140
8		●		HFXI 8-100 Ti L80	HFXI 8-100 Ti L140
9		●		HFXI 9-100 Ti L80	HFXI 9-100 Ti L140
120	6		●	HFXI 6-120 Ti L80	HFXI 6-120 Ti L140
	7	●		HFXI 7-120 Ti L80	HFXI 7-120 Ti L140
	8	●		HFXI 8-120 Ti L80	HFXI 8-120 Ti L140
150	6		●	HFXI 6-150 Ti L80	HFXI 6-150 Ti L140
	7	●		HFXI 7-150 Ti L80	HFXI 7-150 Ti L140
	8	●		HFXI 8-150 Ti L80	HFXI 8-150 Ti L140
200	6		●	HFXI 6-200 Ti L80	HFXI 6-200 Ti L140
	7	●		HFXI 7-200 Ti L80	HFXI 7-200 Ti L140
	8	●		HFXI 8-200 Ti L80	HFXI 8-200 Ti L140

Distributed by

HEXACATH
PIONEER IN BIO ACTIVE COATING

4, Passage Saint Antoine
92508 Rueil-Malmaison Cedex - France
Tel. +33 (0)1 41 39 01 92
Fax. +33 (0)1 41 39 01 93

www.hexacath.com